

Antera 3D® – Comparison Table of Antera Software Versions

| Topography Measurement | Antera 3D | Antera 3D Pro | Antera 3D CS |
|---|-----------|---------------|--------------|
| Wrinkles' overall size (a.u.) | ✓ | ✓ | ✓ |
| Wrinkles' average depth (mm) | ✓ | ✓ | ✓ |
| Wrinkles' average width (mm) | ✓ | ✓ | ✓ |
| Indentation index (a.u.) | ✓ | ✓ | ✓ |
| Wrinkles' maximum depth (mm) | | | ✓ |
| Pores' total volume (mm ³) | | ✓ | ✓ |
| Pores' index (a.u.) | | ✓ | ✓ |
| Pores' count | | ✓ | ✓ |
| Pores' density (/cm ²) | | | ✓ |
| Mean pore volume (mm ³) | | | ✓ |
| Mean pore area (mm ²) | | | ✓ |
| Pores' area (mm ²) | | | ✓ |
| Pores' maximum depth (mm) | | | ✓ |
| Texture – Arithmetical mean roughness Ra (a.u.) | ✓ | ✓ | ✓ |
| Texture - Squared mean roughness Rq (a.u.) | | | ✓ |
| Texture - Elevation span Rt (mm) | | | ✓ |
| Volume of depressions (mm ³) | | ✓ | ✓ |
| Maximum valley depth, Rv (mm) | | | ✓ |
| Volume of elevations (mm ³) | | ✓ | ✓ |
| Maximum peak height, Rp (mm) | | | ✓ |
| Chromophores' Measurement | Antera 3D | Antera 3D Pro | Antera 3D CS |
| Pigmentation average concentration (a.u.) | ✓ | ✓ | ✓ |
| Pigmentation uniformity (a.u.) | ✓ | ✓ | ✓ |
| Pigmentation hyper-concentration (a.u.) | | | ✓ |
| Pigmentation hypo-concentration (a.u.) | | | ✓ |
| Pigmentation hypo-concentration affected area (mm ²) | | | ✓ |
| Redness average concentration (a.u.) | ✓ | ✓ | ✓ |
| Redness uniformity (a.u.) | ✓ | ✓ | ✓ |
| Redness hyper-concentration (a.u.) | | | ✓ |
| Redness hypo-concentration (a.u.) | | | ✓ |
| Redness hypo-concentration affected area (mm ²) | | | ✓ |
| Colour Measurement | Antera 3D | Antera 3D Pro | Antera 3D CS |
| Overall colour variation Δ(L*a*b*) | ✓ | ✓ | ✓ |
| L*a*b* values | | | ✓ |
| L*a*b* variation - Δ L*, Δ a*, Δ b* | | | ✓ |
| Before & After Comparison | Antera 3D | Antera 3D Pro | Antera 3D CS |
| Spot-On | ✓ | ✓ | ✓ |
| Manual matching | ✓ | ✓ | ✓ |
| Manual adjustment of selected area | ✓ | ✓ | ✓ |
| Maximum number of images that can be analysed simultaneously | 10 | 15 | No limit* |
| Ability to save selection | | | ✓ |
| Report | Antera 3D | Antera 3D Pro | Antera 3D CS |
| Create report | ✓ | ✓ | ✓ |
| Save report | ✓ | ✓ | ✓ |
| Print report | ✓ | ✓ | ✓ |
| Export report in HTML and PDF | ✓ | ✓ | ✓ |
| E-mail report | ✓ | ✓ | ✓ |
| Data Export | Antera 3D | Antera 3D Pro | Antera 3D CS |
| Export view in JPEG | ✓ | ✓ | ✓ |
| Export measured parameters | | | ✓ |
| Export raw data (z-elevation, L*a*b* values, pigmentation, redness) | | | ✓ |
| Other Features | Antera 3D | Antera 3D Pro | Antera 3D CS |
| Remote image acquisition | | | ✓ |
| Option to attach camera to a tripod | | | ✓ |

Note: The maximum number of images that can be open depends on the computer memory.